

| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

## **Contents**

| Appendix 2: Evidence of Demonstration Test Results  | 22 |
|---|----|
| Appendix 1: Evidence of Automated Unit Test Results | 15 |
| 8 Approvals   | 14 |
| 7 References  | 13 |
| 6.1 Test procedure                                  | 13 |
| 6 Inspection  | 13 |
| 5 Demonstration Test Results                        | 11 |
| 4 Automated Test Results                            | 2  |
| 3 Responsibilities                                  | 2  |
| 2 Document Scope                                    | 2  |
| 1 Document Purpose                                  | 2  |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

## **1 Document Purpose**

The purpose of this Test Results document is to collect the output of each test defined in OMNI-11-TestProtocols and whether each test has been performed and passed acceptance criteria.

### **2 Document Scope**

This document contains a list of software requirements

Each requirement defined in OMNI-SM-002 Software Requirements Specification is listed in the Software Verification Plan. Requirements are listed in numerical order by requirement number.

The Software Verification Plan includes the following information:

- Software requirement identification number
- Text of the requirement
- Method(s) for verification: Inspection (I), Analysis (A), Demonstration (D), Test (T)

## 3 Responsibilities

| Function            | Responsible For  |
|---------------------|--|
| Software Developers | <ul> <li>Develops and executes verification protocols</li> <li>Writes protocol-specific reports</li> </ul> |
| Lead Developer      | Writes the verification summary report   |

### **4 Automated Test Results**

| Req# | Associated Unit Tests  | Test Result | Initials of<br>Tester and<br>Date Tested |
|------|--|-------------|--|
| 4.3  | SummitServerTest.Startup_CreatesWindow confirms that the microservice creates a window allowing for graceful shutdown. | Passed      | AO<br>11/4/22                            |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| 5.1.1 | SummitServerTest.Startup_PrintsToConsole confirms that the console is written with the current version number and supported devices.  | Passed | AO<br>11/4/22 |
|-------|---|--------|---------------|
|       | Info_Service_InspectRepository confirms that an info service has been created.  | Passed | AO<br>11/4/22 |
|       | ConnectedBridges_WithBridges_ReturnsBridges confirms that a bridge service has been created.  | Passed | AO<br>11/4/22 |
|       | DeviceStatus_WithSuccessfulRead_ReturnsBatteryStatus confirms that a device service has been created.   | Passed | AO<br>11/4/22 |
| 5.1.2 | DeviceStatus_WithSuccessfulRead_ReturnsBatteryStatus confirms that a connected client can request a function to be called on a device instance. This also tests that the function call is forwarded on to the Summit DLL, which after completion of execution returns the response to the client. | Passed | AO<br>11/4/22 |
|       | <b>DeviceStatus_DeviceDoesNotExist_ReturnsError</b> confirms that if the device does not exist, the request is denied.  | Passed | AO<br>11/4/22 |
|       | DescribeBridge_DoesNotExist_ReturnsError confirms that if the bridge does not exist, the request will be denied.  | Passed | AO<br>11/4/22 |
|       | DescribeBridge_NoErrorFromAPI_ReturnsDetails confirms that a connected client can request a function to be called on a bridge instance. This also tests that the function call is forwarded on to the Summit DLL, which after completion of execution returns the response to the client.         | Passed | AO<br>11/4/22 |
| 5.1.3 | SummitServerTest.Startup_CreatesWindow confirms that the microservice creates a window that is shut down by asserting that the process has exited.  | Passed | AO<br>11/4/22 |

| Document is for reference only – not for as-is submission to FDA | Page 3 of 26 |  |
|--|--------------|--|
| ·  | "            |  |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

|       | Shutdown_DisposesManager confirms that device and bridge connections are disposed when the Server is shut down.  | Passed | AO<br>11/4/22 |
|-------|--|--------|---------------|
|       | Startup_CreatesWindow_qToShutdown confirms that the server is shut down when 'q' is entered into the command console.  | Passed | AO<br>11/4/22 |
| 5.2.1 | Info_Service_VersionNumber confirms that the microservice shall return the major, minor, and patch version numbers.  | Passed | AO<br>11/4/22 |
|       | Info_Service_SuppotedDevices confirms that the microservice shall return the list of supported implantable neurostimulators.   | Passed | AO<br>11/4/22 |
|       | Info_Service_InspectRepository confirms that the microservice shall return all cached device and bridge connections.   | Passed | AO<br>11/4/22 |
| 5.2.2 | ListBridges_WithDuplicateBridges_ReturnsUnique, ListBridges_WithQuery_ReturnsFiltered, and ListBridges_WithNoBridges_ReturnsNoBridges confirm that listBridges will return a list of known bridges. Discovery of known bridges is Summit API functionality that is handled through same interface and does not need testing. | Passed | AO<br>11/4/22 |
|       | ConnectedBridges_WithBridges_ReturnsBridges, ConnectedBridges_WithQuery_ReturnsFiltered, and ConnectedBridges_WithNoBridges_ReturnsNoBridges confirm that the microservice shall return a list of bridges currently connected to by the microservice.  | Passed | AO<br>11/4/22 |
|       | DescribeBridge_ErrorFromAPI_ReturnsErrorEmptyDetails and DescribeBridge_NoErrorFromAPI_ReturnsDetails confirm that the microservice shall return telemetry information for a specified bridge  | Passed | AO<br>11/4/22 |
| 5.2.3 | ConnectToBridge_CachedConnectionIsDisposed _ReturnSuccess confirms that when the Summit is in the repo and disposed that the connection call still proceeds.   | Passed | AO<br>11/4/22 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

|       | ConnectToBridge_CachedConnectionNotDisposed _ReturnsSuccess confirms that when the Summit is in the repo but not disposed that the connection call is aborted but success is still returned. | Passed | AO<br>11/4/22 |
|-------|--|--------|---------------|
|       | ConnectToBridge_NoCachedConnectionConnectFails _ReturnsError confirms that Summit connection failures are returned to the application.   | Passed | AO<br>11/4/22 |
|       | ConnectToBridge_NoCachedConnection_ReturnsSuccess _ReturnsSuccess confirms that Summit connection success are returned to the application.   | Passed | AO<br>11/4/22 |
|       | ConnectToBridge_NoDiscoveredBridgeRequest _ReturnsError confirms that an error is returned if an invalid bridge is requested to connect to.  | Passed | AO<br>11/4/22 |
|       | DisconnectFromBridge_DisposesSummit confirms that the DisposeSummit function is called using the mock verify function.   | Passed | AO<br>11/4/22 |
| 5.2.4 | ConfigureBeep_ValidRequestParameter_ReturnsSuccess confirms that the microservice can configure the telemetry beep using API commands.   | Passed | AO<br>11/4/22 |
|       | ConfigureBeep_InvalidRequest_ReturnsError confirms that the microservice returns an error when an invalid parameter is provided to configure the telemetry beep.                             | Passed | AO<br>11/4/22 |
|       | ConfigureBeep_ErrorEncountered_ReturnsFailure confirms that the microservice returns an error when the microservice receives an error from the Summit API.                                   | Passed | AO<br>11/4/22 |
| 5.2.5 | ConnectionStatusStream_NoCachedConnection _ReturnsEmpty confirms that an empty (ended) stream is returned if there is no active connection to monitor.                                       | Passed | AO<br>11/4/22 |

| Document is for reference only – not for as-is submission to FDA Page 5 of 26 | as-is submission to FDA Page 5 of 26 |
|---|--------------------------------------|
|---|--------------------------------------|



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| ConnectionStatusStream_AlreadyConnected_ReturnsEmpty confirms that an empty (ended) stream is returned if there is no active connection to monitor.   | Passed | AO<br>11/4/22 |
|---|--------|---------------|
| ConnectionStatusStream_EnableStreamRequest _ReturnsStream confirms that a stream is returned if a request is made while there is a connection and there is no preexisting connection stream.  | Passed | AO<br>11/4/22 |
| ConnectionStatusStream_ConnectedDisableStream _EndsStream confirms that a stream can be ended by the client requesting to disable the stream.   | Passed | AO<br>11/4/22 |
| ConnectionMonitoring_DisconnectionEventReceivedSummitNotDis  posed_Notify confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a disconnect occurs when the Summit is not disposed. | Passed | AO<br>11/4/22 |
| ConnectionMonitoring_DisconnectionEventReceivedSummitDispose d_Notify confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a disconnect occurs when the Summit is disposed.         | Passed | AO<br>11/4/22 |
| ConnectionMonitoring_DisconnectionEventReceived_Reconnection CTMAttemptsExceeded confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a CTM reconnection attempt fails.             | Passed | AO<br>11/4/22 |
| ConnectionMonitoring_DisconnectionEventReceived_ReconnectionI  NSAttemptsExceeded confirms that the microservice provides an interface for streaming bridge related status that notifies an application when an INS reconnection attempt fails.           | Passed | AO<br>11/4/22 |
| ConnectionMonitoring_DisconnectionEventReceived_Reconnection Succeeded confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a reconnection succeeds.                                | Passed | AO<br>11/4/22 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| ListDevices_WithNoDevices_ReturnsNoDevices confirms that if no devices are discovered by CTM that no devices are returned to the client.   | Passed | AO<br>11/4/22 |
|--|--------|---------------|
| ListDevices_WithSameDeviceConnectedAlready _ReturnsMedtronicError confirms that if a device is already connected to, that the expected Medtronic error code is returned to the client. | Passed | AO<br>11/4/22 |
| ListDevices_WithDeviceSuccess_ReturnsDeviceList confirms that if a device is discovered that it is returned to the client.   | Passed | AO<br>11/4/22 |
| <b>DeviceStatus_WithSuccessfulRead_ReturnsBatteryStatus</b> confirms that on a successful device status call, that battery status information is returned to the client.               | Passed | AO<br>11/4/22 |
| <b>DeviceStatus_DeviceDoesNotExist_ReturnsError</b> confirms that if the device does not exist that an error is returned to the client.  | Passed | AO<br>11/4/22 |
| DeviceStatus_WithUnsuccessfulRead _ReturnsNullBatteryLevels confirms that if an error is encountered during the device status request that an error is returned to the client.         | Passed | AO<br>11/4/22 |
| LeadIntegrityTest_NoError_ReturnsLeadList confirms that the lead integrity test returns impedances if successful.  | Passed | AO<br>11/4/22 |
| <b>LeadIntegrityTest_SameLeadValues_ReturnsZero</b> confirms that the lead integrity test returns 0 if a selected lead pair contains two of the same electrode.                        | Passed | AO<br>11/4/22 |
| LeadIntegrityTest_OutofRangeLeadValues_ReturnsError confirms that the lead integrity test returns out of range error if invalid lead indexes are provided.                             | Passed | AO<br>11/4/22 |
|  |        |               |

| Document is for reference only – not for as-is submission to FDA | Page 7 of 26 |
|--|--------------|
| Document is for reference only – not for as-is submission to FDA | Page 7 of 26 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| 5.2.7 | ConnectToDevice_RequestedDeviceNotCached _ReturnsError confirms that when the requested device is not in the device repo, that the Summit connection is aborted are returned to the application.                        | Passed | AO<br>11/4/22 |
|-------|---|--------|---------------|
|       | ConnectToDevice_NoDeviceSpecified_NullConnectSuccess confirms that if no device is specified in the request, that a null device parameter is provided to the Summit StartInsSession command to initiate a null-connect. | Passed | AO<br>11/4/22 |
|       | ConnectToDevice_CachedDeviceConnect  _ReturnsConnectStatus confirms that the Summit connection StartInsSession command status is returned to the application  | Passed | AO<br>11/4/22 |
|       | ConnectToDevice_NoConnectedBridge _ReturnsError confirms that an error is returned if GetConnectionByName returns null.   | Passed | AO<br>11/4/22 |
|       | DisconnectFromDevice_DisposesSummit confirms that the DisposeSummit function is called using the mock verify function.  | Passed | AO<br>11/4/22 |
| 5.2.8 | ConfigureSense_NullRequestParameter_ReturnsFailure confirms that the microservice returns an error if there is a null parameter provided for part of the sensing configuration.   | Passed | AO<br>11/4/22 |
|       | ConfigureSense_PowerBandEnablesCountIncorrect _ReturnsFailure confirms that an error is returned if the power channel list is too long.   | Passed | AO<br>11/4/22 |
|       | ConfigureSense_NullRequestParametersObject _ReturnsFailure confirms that if an entire sensing configuration object (FFT) is null that an error is returned.   | Passed | AO<br>11/4/22 |
|       | ConfigureSense_WriteSensingFftSettingsReject _ReturnsFailure that if a command is rejected as part of the sense configuration transaction that the expected summit API reject is returned to the client.                | Passed | AO<br>11/4/22 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

|       | ConfigureSense_DisabledTimeDomainStream _SetsSampleRateToDisabled confirms that if a time domain stream is indicated to be disabled that the correct enumeration for disabled channels is transmitted to the Summit API. | Passed | AO<br>11/4/22  |
|-------|--|--------|----------------|
|       | ConfigureSense_ValidConfiguration_Success confirms that the microservice returns success when no error is encountered in the sense configuration transaction.  | Passed | AO<br>11/4/22  |
| 5.2.9 | ConfigureSense_NullRequestParameter_ReturnsFailure confirms that the microservice returns an error if there is a null parameter provided for part of the sensing configuration.  | Passed | AO<br>11/11/22 |
|       | ConfigureSense_PowerBandEnablesCountIncorrect _ReturnsFailure confirms that an error is returned if the power channel list is too long.  | Passed | AO<br>11/11/22 |
|       | ConfigureSense_NullRequestParametersObject _ReturnsFailure confirms that if an entire sensing configuration object (FFT) is null that an error is returned.  | Passed | AO<br>11/11/22 |
|       | ConfigureSense_WriteSensingFftSettingsReject _ReturnsFailure that if a command is rejected as part of the sense configuration transaction that the expected summit API reject is returned to the client.                 | Passed | AO<br>11/11/22 |
|       | ConfigureSense_DisabledTimeDomainStream _SetsSampleRateToDisabled confirms that if a time domain stream is indicated to be disabled that the correct enumeration for disabled channels is transmitted to the Summit API. | Passed | AO<br>11/11/22 |
|       | ConfigureSense_ValidConfiguration_Success confirms that the microservice returns success when no error is encountered in the sense configuration transaction.  | Passed | AO<br>11/11/22 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| 5.3.1 | ConnectionMonitoring_DisconnectionEventReceivedSummitNotDis  posed_Notify confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a disconnect occurs when the Summit is not disposed. | Passed | AO<br>11/4/22 |
|-------|---|--------|---------------|
|       | ConnectionMonitoring_DisconnectionEventReceivedSummitDispose d_Notify confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a disconnect occurs when the Summit is disposed.         | Passed | AO<br>11/4/22 |
|       | ConnectionMonitoring_DisconnectionEventReceived_Reconnection CTMAttemptsExceeded confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a CTM reconnection attempt fails.             | Passed | AO<br>11/4/22 |
|       | ConnectionMonitoring_DisconnectionEventReceived_ReconnectionI  NSAttemptsExceeded confirms that the microservice provides an interface for streaming bridge related status that notifies an application when an INS reconnection attempt fails.           | Passed | AO<br>11/4/22 |
|       | ConnectionMonitoring_DisconnectionEventReceived_Reconnection Succeeded confirms that the microservice provides an interface for streaming bridge related status that notifies an application when a reconnection suceeds.                                 | Passed | AO<br>11/4/22 |
| 6.1.1 | MultipleConnections_TwoClientOneDeviceConnectionRequestsSuccess and MultipleConnections_TwoClientTwoDevicesConnectionRequestsSuccess confirm that multiple clients can successfully communicate with microservice-connected devices.                      | Passed | AO<br>11/4/22 |
| 6.1.2 | Create <datatype>StreamStream_AlreadyConnected _ReturnsEmpty confirms that an error is returned to the client if the specific data type stream is already active.</datatype>  | Passed | AO<br>11/4/22 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| 6.2.1 | MultipleConnections_TwoClientTwoDevices                                | Passed | AO      |
|-------|--|--------|---------|
|       | _ConnectionRequestsSuccess and   |        | 11/4/22 |
|       | MultipleConnections_OneClientTwoDevices                                |        |         |
|       | _ConnectionRequestsSuccess confirm that the microservice allows        |        |         |
|       | for two simultaneous connections to Summit devices.                    |        |         |
| 8.4.1 | ListBridges_WithDuplicateBridges_ReturnsUnique() confirms that         | Passed | AO      |
|       | without any configuration, localhost is used to host the microservice. |        | 11/4/22 |
| 10.4  | Generation of signed interface control document (ICD) demonstrates     | Passed | AO      |
|       | completion of this requirement.  |        | 11/4/22 |

## **5 Demonstration Test Results**

| Req# | Associated Demonstration Protocol   | Test Result   | Initials of<br>Tester<br>and Date<br>Tested |
|------|---|---|---|
| ex.1 | Example 1 – take screenshot of system running   | Passed See screenshot in Appendix 1 CTM: NKWXXXXX INS: NPCXXXXX     | DAB<br>2022-11-02                           |
| 4.2  | <ol> <li>A tester will run the microservice<br/>connected to a Medtronic Summit System<br/>and verify that they are connected.</li> </ol>   | Passed See screenshot in Appendix 1 CTM: NKW021653N INS: NPC700365H | AO<br>2022-11-03                            |
| 7.1  | <ol> <li>A tester will download the microservice executable and run it on a Windows 10 64 bit PC.</li> <li>The tester will take a screenshot of the console window running on the PC and upload to the test report document.</li> </ol> | Passed See screenshot in Appendix 1 CTM: NKW021653N INS: NPC700365H | AO<br>2022-11-03                            |

| Document is for reference only – not for as-is submission to FDA | Page 11 of 26  |
|--|----------------|
| boddinent is for reference only mortal as is submission to 15%   | 1 486 11 01 20 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| 7.2  | <ol> <li>A tester will run the microservice using gRPC version 1.0.0 or greater and verify that after following all of the setup steps above that the application runs successfully and connects to the CTM and Summit.</li> <li>The tester will document in the test report what version of gRPC the client application utilized when the demonstration was executed. This version has to be greater than 1.0 for the test to be considered a success.</li> </ol> | Passed<br>gRPC version: 2.36.0<br>See screenshot in Appendix 1<br>CTM: NKW021653N<br>INS: NPC700365H | AO<br>2022-11-03 |
|------|--|--|------------------|
| 7.3  | <ol> <li>A tester will verify that messages are serialized using protobuf version 3. The tester can verify this by checking the version of protobuf</li> <li>The tester will document in the test report what version of protobuf the client application utilized when the demonstration was executed. This version has to be greater than 3.0 for the test to be considered a success.</li> </ol>   | Passed Protobuf version: 3.15.0 CTM: NKW021653N INS: NPC700365H                                      | AO<br>2022-11-03 |
| 9.1  | A tester will verify that the provided installer creates the OmniSummitMicroservice binary. Running the binary verifies that the Medtronic Summit DLLs are in the correct location.  | Passed<br>CTM: NKW021653N<br>INS: NPC700365H   | AO<br>2022-11-03 |
| 10.1 | A tester will look over the version history and verify that it follows proper semantic versioning  | Passed<br>CTM: NKW021653N<br>INS: NPC700365H   | AO<br>2022-11-03 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

| 10.2 | 1) A tester will verify that the installer's   | Passed          | AO         |
|------|--|-----------------|------------|
|      | version reflects the version of the            | CTM: NKW021653N | 2022-11-03 |
|      | microservice                                   | INS: NPC700365H |            |
| 10.3 | 1) A tester will verify that there is a github | Passed          | AO         |
|      | action in place to update and automatically    | CTM: NKW021653N | 2022-11-03 |
|      | deploy documentation for each new              | INS: NPC700365H |            |
|      | release  |                 |            |

# **6 Inspection**

### **6.1 Test procedure**

There are no requirements that are tested via inspection.

## **7 References**

| Document Identifier | Title   |
|---------------------|---|
| OMNI-SM-002         | Software Requirements Specification, OMNI Summit Microservice |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 2    |

# 8 Approvals

| Approver Role  | Signature and Date  |
|----------------|---|
| Project Lead   | Down Borton<br>11/630/1202246E1:54 PM EST                     |
| Quality        | Docusigned by: Timothy Denison 12/0393500325B4FA.6:15 AM EST  |
| Lead Developer | Docusigned by:  Jeffrey Herron  11/3097210722542511:08 AM PST |
| Author         | Docusigned by:  Are Oulsner  11/30750927648/9:58 PM EST       |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

# **Appendix 1: Evidence of Automated Unit Test Results**

Below is the direct output of the console for all unit tests.

| Req# | Evidence  |
|------|---|
|      | Test run for D:\a\OmniSummitDeviceService\OmniSummitDeviceService\OmniSummitDeviceService.UnitTests\bin\Debug\n et472\OmniSummitDeviceService.UnitTests.dll (.NETFramework,Version=v4.7.2) Microsoft (R) Test Execution Command Line Tool Version 16.11.0 Copyright (c) Microsoft Corporation. All rights reserved.   |
|      | Starting test execution, please wait A total of 1 test files matched the specified pattern. NUnit Adapter 4.0.0.0: Test execution started Running all tests in D:\a\OmniSummitDeviceService\OmniSummitDeviceService\OmniSummitDeviceService.UnitTests\bin\Debug\n et472\OmniSummitDeviceService.UnitTests.dll NUnit3TestExecutor discovered 92 of 92 NUnit test cases using Current Discovery mode, Non-Explicit run 721360 |

| Document is for reference only – not for as-is submission to FDA | Page 15 of 26 |
|--|---------------|
| Document is for reference only – not for as-is submission to FDA | Page 15 01 26 |



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

```
Passed Startup CreatesWindow [4 s]
 Passed Startup CreatesWindow qToShutdown [5 s]
OpenMindServer started on port localhost:50051
Version: 1.10.2.0
Supported devices: Medtronic Summit RC+S
Press 'q' key to stop server...
OpenMindServer started on port localhost:50051
Version: 1.10.2.0
Supported devices: Medtronic Summit RC+S
Press 'q' key to stop server...
 Passed Startup PrintsToConsole [1 s]
 Passed ConfigureBeep ErrorEncountered ReturnsError [5 s]
{ "name": "//summit/bridge/foo/device/foo", "error": { "errorType": "API_ERROR", "message":
"NoError", "txTime": "2022-11-17T18:57:12.573098800Z", "rxTime": "2022-11-17T18:57:12.573098800Z",
"linkStatus": "AA==" } }
 Passed ConfigureBeep InvalidRequest ReturnsError [83 ms]
 Passed ConfigureBeep ValidRequestParameter ReturnsSuccess [29 ms]
 Passed ConfigureSense_DisabledTimeDomainStream_SetsSampleRateToDisabled [338 ms]
  Passed ConfigureSense NullRequestParameter ReturnsFailure [20 ms]
  Passed ConfigureSense NullRequestParametersObject ReturnsFailure [25 ms]
```



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

```
Passed ConfigureSense PowerBandConfigurationCountIncorrect ReturnsFailure [18 ms]
Passed ConfigureSense PowerBandEnablesCountIncorrect ReturnsFailure [51 ms]
Passed ConfigureSense ValidConfiguration Success [639 ms]
Passed ConfigureSense_WriteSensingFftSettingsReject_ReturnsFailure [25 ms]
Passed ConnectedBridges_WithBridges_ReturnsBridges [25 ms]
Passed ConnectedBridges_WithNoBridges_ReturnsNoBridges [22 ms]
Passed ConnectedBridges_WithQuery_ReturnsFiltered [14 ms]
Passed ConnecteToDevice NoDeviceSpecified NullConnectSuccess [67 ms]
Passed ConnectionMonitoring DisconnectionEventReceived ReconnectionCTMAttemptsExceeded [20 s]
Passed ConnectionMonitoring DisconnectionEventReceived ReconnectionINSAttemptsExceeded [15 s]
Passed ConnectionMonitoring DisconnectionEventReceived ReconnectionSucceeded [5 s]
Passed ConnectionMonitoring DisconnectionEventReceivedSummitDisposed Notify [137 ms]
Passed ConnectionMonitoring DisconnectionEventReceivedSummitNotDisposed Notify [162 ms]
Passed ConnectionStatusStream_AlreadyConnected_ReturnsEmpty [119 ms]
Passed ConnectionStatusStream ConnectedDisableStream EndsStream [132 ms]
Passed ConnectionStatusStream EnableStreamRequest ReturnsStream [132 ms]
Passed ConnectionStatusStream NoCachedConnection ReturnsEmpty [28 ms]
Passed ConnectToBridge_CachedConnectionIsDisposed_ReturnsSuccess [23 ms]
Passed ConnectToBridge CachedConnectionNotDisposed ReturnsSuccess [24 ms]
Passed ConnectToBridge NoCachedConnectionConnectFails ReturnsError [19 ms]
Passed ConnectToBridge_NoCachedConnectionRepoSuccess_ReturnsSuccess [28 ms]
Passed ConnectToBridge NoDiscoveredBridgeRequest ReturnsError [19 ms]
Passed ConnectToDevice CachedDeviceConnect ReturnsConnectStatus [21 ms]
```



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

```
Passed ConnectToDevice NoConnectedBridge ReturnsError [21 ms]
Passed ConnectToDevice RequestedDeviceNotCached ReturnsError [20 ms]
Passed CreateAdaptiveDataStream AlreadyConnected ReturnsEmpty [119 ms]
Passed CreateAdaptiveDataStream ConnectedDisableStream EndsStream [124 ms]
Passed CreateAdaptiveDataStream_EnableStreamRequestData_ReturnsStream [173 ms]
Passed CreateAdaptiveDataStream NoCachedConnection ReturnsEmpty [26 ms]
Passed CreateBandPowerDataStream_AlreadyConnected_ReturnsEmpty [123 ms]
Passed CreateBandPowerDataStream ConnectedDisableStream EndsStream [131 ms]
Passed CreateBandPowerDataStream EnableStreamRequestData ReturnsStream [145 ms]
Passed CreateBandPowerDataStream NoCachedConnection ReturnsEmpty [27 ms]
Passed CreateEchoDataStream AlreadyConnected ReturnsEmpty [119 ms]
Passed CreateEchoDataStream ConnectedDisableStream EndsStream [129 ms]
Passed CreateEchoDataStream EnableStreamRequestData ReturnsStream [138 ms]
Passed CreateEchoDataStream NoCachedConnection ReturnsEmpty [19 ms]
Passed CreateFourierDataStream AlreadyConnected ReturnsEmpty [123 ms]
Passed CreateFourierDataStream ConnectedDisableStream EndsStream [128 ms]
Passed CreateFourierDataStream_EnableStreamRequestData_ReturnsStream [147 ms]
Passed CreateFourierDataStream_NoCachedConnection_ReturnsEmpty [26 ms]
Passed CreateInertialDataStream AlreadyConnected ReturnsEmpty [117 ms]
Passed CreateInertialDataStream ConnectedDisableStream EndsStream [127 ms]
Passed CreateInertialDataStream_EnableStreamRequestData ReturnsStream [180 ms]
Passed CreateInertialDataStream NoCachedConnection ReturnsEmpty [22 ms]
Passed CreateLoopDataStream AlreadyConnected ReturnsEmpty [121 ms]
```



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

```
Passed CreateLoopDataStream ConnectedDisableStream EndsStream [132 ms]
Passed CreateLoopDataStream EnableStreamRequestData ReturnsStream [125 ms]
Passed CreateLoopDataStream NoCachedConnection ReturnsEmpty [24 ms]
Passed CreateTimeDomainDataStream AlreadyConnected ReturnsEmpty [125 ms]
Passed CreateTimeDomainDataStream_ConnectedDisableStream_EndsStream [130 ms]
Passed CreateTimeDomainDataStream_EnableStreamRequest_ReturnsStream [143 ms]
Passed CreateTimeDomainDataStream_NoCachedConnection_ReturnsEmpty [21 ms]
Passed DescribeBridge DoesNotExist ReturnsError [22 ms]
Passed DescribeBridge ErrorFromAPI ReturnsErrorEmptyDetails [24 ms]
Passed DescribeBridge NoErrorFromAPI ReturnsDetails [19 ms]
Passed DeviceStatus DeviceDoesNotExist ReturnsError [23 ms]
Passed DeviceStatus WithSuccessfulRead ReturnsBatteryStatus [24 ms]
Passed DeviceStatus WithUnsuccessfulRead ReturnsNullBatteryLevels [22 ms]
Passed DisconnectFromBridge DisposesSummit [23 ms]
Passed DisconnectFromDevice DisposesSummit [22 ms]
Passed Info Service InspectRepository [21 ms]
Passed Info Service SuppotedDevices [22 ms]
Passed Info_Service_VersionNumber [18 ms]
Passed LeadIntegrityTest NoError ReturnsLeadList [27 ms]
Passed LeadIntegrityTest OutofRangeLeadValues ReturnsError [20 ms]
Passed LeadIntegrityTest SameLeadValues ReturnsZero [24 ms]
Passed ListBridges WithDuplicateBridges ReturnsUnique [19 ms]
Passed ListBridges WithNoBridges ReturnsNoBridges [21 ms]
```



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

```
Passed ListBridges WithQuery ReturnsFiltered [18 ms]
  Passed ListDevices WithDeviceSuccess ReturnsDeviceList [30 ms]
  Passed ListDevices WithNoDevices ReturnsNoDevices [19 ms]
  Passed ListDevices WithSameDeviceConnectedAlready ReturnsMedtronicError [21 ms]
  Passed MultipleConnections_OneClientTwoDevices_ConnectionRequestsSuccess [21 ms]
NUnit Adapter 4.0.0.0: Test execution complete
  Passed MultipleConnections_TwoClientOneDevice_ConnectionRequestsSuccess [23 ms]
  Passed MultipleConnections TwoClientTwoDevices ConnectionRequestsSuccess [20 ms]
  Passed Shutdown DisposesManager [1 ms]
  Passed StreamDisable_ErrorEncountered ReturnsFailure [27 ms]
  Passed StreamDisable InvalidRequestParameter ReturnsUnknown [17 ms]
  Passed StreamDisable ValidRequestParameter ReturnsSuccess [23 ms]
  Passed StreamEnable_ErrorEncountered_ReturnsFailure [20 ms]
 Passed StreamEnable InvalidRequestParameter ReturnsUnknown [20 ms]
  Passed StreamEnable ValidRequestParameter ReturnsSuccess [19 ms]
Test Run Successful.
Total tests: 92
     Passed: 92
 Total time: 1.3951 Minutes
     1>Done Building Project
"D:\a\OmniSummitDeviceService\OmniSummitDeviceService\OmniSummitDeviceService.sln" (VSTest
target(s)).
```



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

3>Done Building Project

"D:\a\OmniSummitDeviceService\OmniSummitDeviceService\OmniSummitDeviceService.UnitTests\OmniSummit DeviceService.UnitTests.csproj" (VSTest target(s)).

Build succeeded.

- 0 Warning(s)
- 0 Error(s)

Time Elapsed 00:01:28.04




| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

# **Appendix 2: Evidence of Demonstration Test Results**

| Req# |
|------|
|------|



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |




Document is for reference only – not for as-is submission to FDA

Page 23 of 26




| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |





| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |



Document is for reference only – not for as-is submission to FDA

Page 25 of 26



| Doc. Number | Rev. |
|-------------|------|
| OMNI-SM-012 | 1    |

Document is for reference only – not for as-is submission to FDA